

Sortilin-related receptor

Catalog No: #AP87935



Package Size: #AP87935-1 50ug #AP87935-2 100ug #AP87935-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sortilin-related receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Low-density lipoprotein receptor relative with 11 ligand-binding repeats,LDLR relative with 11 ligand-binding repeats,SorLA-1,Sorting protein-related receptor containing LDLR class A repeats,SorLA,C11orf32
Accession No.	Uniprot ID:Q92673
Uniprot	Q92673
GeneID	6653;
SDS-PAGE MW	15.84
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:227-370aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SORL1

Note: This product is for in vitro research use only