

## D-3-phosphoglycerate dehydrogenase

Catalog No: #AP87936



Package Size: #AP87936-1 50ug #AP87936-2 100ug #AP87936-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | D-3-phosphoglycerate dehydrogenase                     |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Bovine   |
| Other Names        | 3-PGDH,2-oxoglutarate reductase,Malate dehydrogenase   |
| Accession No.      | Uniprot ID:Q5EAD2                                      |
| Uniprot            | Q5EAD2   |
| GeneID             | 505103;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PHGDH

Note: This product is for in vitro research use only