

# SPARC

Catalog No: #AP87938



Package Size: #AP87938-1 50ug #AP87938-2 100ug #AP87938-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	SPARC
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Basement-membrane protein 40,BM-40,Osteonectin,ON,Secreted protein acidic and rich in cysteine
Accession No.	Uniprot ID:P36377
Uniprot	P36377
GeneID	386571;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SPARC

Note: This product is for in vitro research use only