

SPARC

Catalog No: #AP87940



Package Size: #AP87940-1 50ug #AP87940-2 100ug #AP87940-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	SPARC
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Basement-membrane protein 40,BM-40,Osteonectin,ON,Secreted protein acidic and rich in cysteine
Accession No.	Uniprot ID:P20112
Uniprot	P20112
GeneID	595124;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SPARC

Note: This product is for in vitro research use only