

Histidine--tRNA ligase, cytoplasmic

Catalog No: #AP87942



Package Size: #AP87942-1 50ug #AP87942-2 100ug #AP87942-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Histidine--tRNA ligase, cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Histidyl-tRNA synthetase,HisRS
Accession No.	Uniprot ID:Q61035
Uniprot	Q61035
GeneID	15115;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hars

Note: This product is for in vitro research use only