

## Serine/threonine-protein kinase STK11

Catalog No: #AP87949



Package Size: #AP87949-1 50ug #AP87949-2 100ug #AP87949-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine/threonine-protein kinase STK11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Liver kinase B1,LKB1,Renal carcinoma antigen NY-REN-19,LKB1,PJS
Accession No.	Uniprot ID:Q15831
Uniprot	Q15831
GeneID	6794;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:STK11

Note: This product is for in vitro research use only