

Signal transducer and activator of transcription 6

Catalog No: #AP87950



Package Size: #AP87950-1 50ug #AP87950-2 100ug #AP87950-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Signal transducer and activator of transcription 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-4 Stat
Accession No.	Uniprot ID:P42226
Uniprot	P42226
GeneID	6778;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: STAT6

Note: This product is for in vitro research use only