

Protein S100-A12

Catalog No: #AP87957



Package Size: #AP87957-1 50ug #AP87957-2 100ug #AP87957-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein S100-A12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Calcium-binding protein in amniotic fluid 1, CAAF1, Calgranulin-C, CAGC, Cornea-associated antigen, CO-AG, Extracellular newly identified RAGE-binding protein, EN-RAGE, RAGE-binding protein, S100 calcium-binding protein A12, CAAF1
Accession No.	Uniprot ID:P79105
Uniprot	P79105
GeneID	282467;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: S100A12

Note: This product is for in vitro research use only