

## Long-chain fatty acid transport protein 1

Catalog No: #AP87958



Package Size: #AP87958-1 50ug #AP87958-2 100ug #AP87958-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Long-chain fatty acid transport protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FATP-1,Solute carrier family 27 member 1,ACSVL5,FATP1
Accession No.	Uniprot ID:Q6PCB7
Uniprot	Q6PCB7
GeneID	376497;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SLC27A1

Note: This product is for in vitro research use only