

Solute carrier family 40 member 1

Catalog No: #AP87959



Package Size: #AP87959-1 50ug #AP87959-2 100ug #AP87959-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Solute carrier family 40 member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ferroportin-1, Iron-regulated transporter 1, Metal transporter protein 1, MTP1, Fpn1, Ireg1, Slc11a3, Slc39a1
Accession No.	Uniprot ID:Q9JHI9
Uniprot	Q9JHI9
GeneID	53945;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Slc40a1

Note: This product is for in vitro research use only