

MAGB3 Antibody

Catalog No: #44468

Package Size: #44468-1 50ul #44468-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MAGB3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	MAGB3 Antibody detects endogenous levels of total MAGB3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MAGB3
Target Name	MAGB3
Other Names	MAGEB3, CT3.5, MAGE-B3 antigen, Melanoma-associated antigen B3, Melanoma antigen family B, 3
Accession No.	Swiss-Prot#: O15480NCBI Gene ID:4114
Uniprot	O15480
GeneID	4114;
Calculated MW	39 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only