

S-phase kinase-associated protein 2

Catalog No: #AP87961



Package Size: #AP87961-1 50ug #AP87961-2 100ug #AP87961-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	S-phase kinase-associated protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cyclin-A/CDK2-associated protein p45,F-box protein Skp2,F-box/WD-40 protein 1,FWD1
Accession No.	Uniprot ID:Q9Z0Z3
Uniprot	Q9Z0Z3
GeneID	27401;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Skp2

Note: This product is for in vitro research use only