

## Signal transducer and activator of transcription 3

Catalog No: #AP87965



Package Size: #AP87965-1 50ug #AP87965-2 100ug #AP87965-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Signal transducer and activator of transcription 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Acute-phase response factor, Aprf
Accession No.	Uniprot ID:P42227
Uniprot	P42227
GeneID	20848;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Stat3

Note: This product is for in vitro research use only