

Bone sialoprotein 2

Catalog No: #AP87966



Package Size: #AP87966-1 50ug #AP87966-2 100ug #AP87966-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bone sialoprotein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Bone sialoprotein II,BSP II,Cell-binding sialoprotein,Integrin-binding sialoprotein
Accession No.	Uniprot ID:P79780
Uniprot	P79780
GeneID	396071;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IBSP

Note: This product is for in vitro research use only