

Sclerostin domain-containing protein 1

Catalog No: #AP87970



Package Size: #AP87970-1 50ug #AP87970-2 100ug #AP87970-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sclerostin domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ectodermal BMP inhibitor,Ectodin,Sclerostin-like protein,Uterine sensitization-associated gene 1 protein,USAG-1,Sostl,Usag1
Accession No.	Uniprot ID:Q9CQN4
Uniprot	Q9CQN4
GeneID	66042;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sostdc1

Note: This product is for in vitro research use only