

## D-3-phosphoglycerate dehydrogenase

Catalog No: #AP87972



Package Size: #AP87972-1 50ug #AP87972-2 100ug #AP87972-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	D-3-phosphoglycerate dehydrogenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	3-PGDH,2-oxoglutarate reductase,Malate dehydrogenase,PGDH3
Accession No.	Uniprot ID:O43175
Uniprot	O43175
GeneID	26227;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PHGDH

Note: This product is for in vitro research use only