

NAD-dependent deacetylase sirtuin-3, mitochondrial

Catalog No: #AP87973



Package Size: #AP87973-1 50ug #AP87973-2 100ug #AP87973-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NAD-dependent deacetylase sirtuin-3, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	hSIRT3, SIR2-like protein 3, SIR2L3
Accession No.	Uniprot ID:Q9NTG7
Uniprot	Q9NTG7
GeneID	23410;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SIRT3

Note: This product is for in vitro research use only