

Slit homolog 2 protein

Catalog No: #AP87974



Package Size: #AP87974-1 50ug #AP87974-2 100ug #AP87974-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Slit homolog 2 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Slit-2
Accession No.	Uniprot ID:Q9R1B9
Uniprot	Q9R1B9
GeneID	20563;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Slit2

Note: This product is for in vitro research use only