

## Suppressor of cytokine signaling 3

Catalog No: #AP87975



Package Size: #AP87975-1 50ug #AP87975-2 100ug #AP87975-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Suppressor of cytokine signaling 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	SOCS-3,Cytokine-inducible SH2 protein 3,Cish3
Accession No.	Uniprot ID:O88583
Uniprot	O88583
GeneID	89829;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Socs3

Note: This product is for in vitro research use only