

Sphingosine 1-phosphate receptor 4

Catalog No: #AP87976



Package Size: #AP87976-1 50ug #AP87976-2 100ug #AP87976-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sphingosine 1-phosphate receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	S1P receptor 4, Endothelial differentiation G-protein coupled receptor 6, Lysophospholipid receptor C1, Sphingosine 1-phosphate receptor Edg-6, S1P receptor Edg-6, Edg6, Lpc1, S1p4
Accession No.	Uniprot ID: Q9Z0L1
Uniprot	Q9Z0L1
GeneID	13611;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: S1pr4

Note: This product is for in vitro research use only