

Sphingosine 1-phosphate receptor 2

Catalog No: #AP87980



Package Size: #AP87980-1 50ug #AP87980-2 100ug #AP87980-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sphingosine 1-phosphate receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	S1P receptor 2,AGR16,Endothelial differentiation G-protein coupled receptor 5,G-protein coupled receptor H218,Sphingosine 1-phosphate receptor Edg-5,S1P receptor Edg-5,Edg5
Accession No.	Uniprot ID:P47752
Uniprot	P47752
GeneID	29415;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:S1pr2

Note: This product is for in vitro research use only