

Pendrin

Catalog No: #AP87981



Package Size: #AP87981-1 50ug #AP87981-2 100ug #AP87981-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pendrin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Sodium-independent chloride/iodide transporter,Solute carrier family 26 member 4,Pds
Accession No.	Uniprot ID:Q9R154
Uniprot	Q9R154
GeneID	29440;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc26a4

Note: This product is for in vitro research use only