

Beta-galactoside alpha-2,6-sialyltransferase 1

Catalog No: #AP87982



Package Size: #AP87982-1 50ug #AP87982-2 100ug #AP87982-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-galactoside alpha-2,6-sialyltransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alpha 2,6-ST 1,CMP-N-acetylneuraminate-beta-galactosamide-alpha-2,6-sialyltransferase 1,ST6Gal I,ST6GalI,Sialyltransferase 1,Siat1
Accession No.	Uniprot ID:Q64685
Uniprot	Q64685
GeneID	20440;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:St6gal1

Note: This product is for in vitro research use only