

## Sialic acid-binding Ig-like lectin 10

Catalog No: #AP87983



Package Size: #AP87983-1 50ug #AP87983-2 100ug #AP87983-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sialic acid-binding Ig-like lectin 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Siglec-10, Siglec-like protein 2, SLG2, UNQ477/PRO940
Accession No.	Uniprot ID:Q96LC7
Uniprot	Q96LC7
GeneID	89790;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: SIGLEC10

Note: This product is for in vitro research use only