

Smoothened homolog

Catalog No: #AP87984



Package Size: #AP87984-1 50ug #AP87984-2 100ug #AP87984-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Smoothened homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	SMO,SMOH
Accession No.	Uniprot ID:O42224
Uniprot	O42224
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SMO

Note: This product is for in vitro research use only