

NAD-dependent deacetylase sirtuin-2

Catalog No: #AP87986



Package Size: #AP87986-1 50ug #AP87986-2 100ug #AP87986-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NAD-dependent deacetylase sirtuin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	SIR2-like protein 2, Sir2I2
Accession No.	Uniprot ID:Q5RJQ4
Uniprot	Q5RJQ4
GeneID	361532;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Sirt2

Note: This product is for in vitro research use only