

Long-chain fatty acid transport protein 1

Catalog No: #AP87988



Package Size: #AP87988-1 50ug #AP87988-2 100ug #AP87988-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Long-chain fatty acid transport protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	FATP-1,Solute carrier family 27 member 1,FATP1
Accession No.	Uniprot ID:Q3ZKN0
Uniprot	Q3ZKN0
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLC27A1

Note: This product is for in vitro research use only