

Tyrosine-protein phosphatase non-receptor type 1

Catalog No: #AP87994



Package Size: #AP87994-1 50ug #AP87994-2 100ug #AP87994-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tyrosine-protein phosphatase non-receptor type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Protein-tyrosine phosphatase 1B,PTP-1B,Protein-tyrosine phosphatase HA2,PTP-HA2
Accession No.	Uniprot ID:P35821
Uniprot	P35821
GeneID	19246;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ptpn1

Note: This product is for in vitro research use only