

TATA-box-binding protein

Catalog No: #AP87999



Package Size: #AP87999-1 50ug #AP87999-2 100ug #AP87999-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TATA-box-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TATA sequence-binding protein,TATA-binding factor,TATA-box factor,Transcription initiation factor TFIID TBP subunit,GTF2D1,TF2D,TFIID
Accession No.	Uniprot ID:P20226
Uniprot	P20226
GeneID	6908;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TBP

Note: This product is for in vitro research use only