

Transferrin receptor protein 1

Catalog No: #AP88000



Package Size: #AP88000-1 50ug #AP88000-2 100ug #AP88000-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transferrin receptor protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TR,T9,p90,CD71
Accession No.	Uniprot ID:P02786
Uniprot	P02786
GeneID	7037;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TFRC

Note: This product is for in vitro research use only