

Taste receptor type 2 member 38

Catalog No: #AP88008



Package Size: #AP88008-1 50ug #AP88008-2 100ug #AP88008-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Taste receptor type 2 member 38
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	T2R38,PTC bitter taste receptor,Taste receptor type 2 member 61,T2R61,PTC
Accession No.	Uniprot ID:P59533
Uniprot	P59533
GeneID	5726;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TAS2R38

Note: This product is for in vitro research use only