

Sodium-dependent noradrenaline transporter

Catalog No: #AP88009



Package Size: #AP88009-1 50ug #AP88009-2 100ug #AP88009-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sodium-dependent noradrenaline transporter
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Norepinephrine transporter,NET,Solute carrier family 6 member 2,NAT1,NET1,SLC6A5
Accession No.	Uniprot ID:P23975
Uniprot	P23975
GeneID	6530;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLC6A2

Note: This product is for in vitro research use only