

Mothers against decapentaplegic homolog 4

Catalog No: #AP88013



Package Size: #AP88013-1 50ug #AP88013-2 100ug #AP88013-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mothers against decapentaplegic homolog 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MAD homolog 4, SMAD family member 4, SMAD 4, Madh4
Accession No.	Uniprot ID:O70437
Uniprot	O70437
GeneID	50554;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Smad4

Note: This product is for in vitro research use only