

Synapsin-1

Catalog No: #AP88015



Package Size: #AP88015-1 50ug #AP88015-2 100ug #AP88015-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Synapsin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Synapsin I, Syn-1
Accession No.	Uniprot ID:O88935
Uniprot	O88935
GeneID	20964;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Syn1

Note: This product is for in vitro research use only