

Zinc transporter ZIP4

Catalog No: #AP88016



Package Size: #AP88016-1 50ug #AP88016-2 100ug #AP88016-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Zinc transporter ZIP4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Solute carrier family 39 member 4,Zrt- and Irt-like protein 4,ZIP-4,ZIP4
Accession No.	Uniprot ID:Q1KZG0
Uniprot	Q1KZG0
GeneID	510951;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLC39A4

Note: This product is for in vitro research use only