

Sphingosine-1-phosphate lyase 1

Catalog No: #AP88017

Package Size: #AP88017-1 50ug #AP88017-2 100ug #AP88017-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sphingosine-1-phosphate lyase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	S1PL, Sphingosine-1-phosphate aldolase, Spl
Accession No.	Uniprot ID:Q8CHN6
Uniprot	Q8CHN6
GeneID	286896;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Sgpl1

Note: This product is for in vitro research use only