

Sphingosine-1-phosphate phosphatase 1

Catalog No: #AP88018

Package Size: #AP88018-1 50ug #AP88018-2 100ug #AP88018-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sphingosine-1-phosphate phosphatase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SPPase1,Sphingosine-1-phosphatase 1
Accession No.	Uniprot ID:Q9BX95
Uniprot	Q9BX95
GeneID	81537;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SGPP1

Note: This product is for in vitro research use only