

Transmembrane protein 126A

Catalog No: #AP88022



Package Size: #AP88022-1 50ug #AP88022-2 100ug #AP88022-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transmembrane protein 126A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:Q9H061
Uniprot	Q9H061
GeneID	84233;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TMEM126A

Note: This product is for in vitro research use only