

Estrogen sulfotransferase, isoform 3

Catalog No: #AP88026



Package Size: #AP88026-1 50ug #AP88026-2 100ug #AP88026-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Estrogen sulfotransferase, isoform 3 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | EST-3,Estrone sulfotransferase,Sulfotransferase, estrogen-preferring |
| Accession No. | Uniprot ID:P49889 |
| Uniprot | P49889 |
| GeneID | 360268; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ste

Note: This product is for in vitro research use only