

Signal transducer and activator of transcription 3

Catalog No: #AP88027



Package Size: #AP88027-1 50ug #AP88027-2 100ug #AP88027-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Signal transducer and activator of transcription 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	RCJMB04_38I20
Accession No.	Uniprot ID:Q6DV79
Uniprot	Q6DV79
GeneID	420027;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: STAT3

Note: This product is for in vitro research use only