

## Transmembrane protein 132D

Catalog No: #AP88028

Package Size: #AP88028-1 50ug #AP88028-2 100ug #AP88028-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Transmembrane protein 132D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Mature oligodendrocytes transmembrane protein,Mature OL transmembrane protein,Molt
Accession No.	Uniprot ID:Q76HP3
Uniprot	Q76HP3
GeneID	243274;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Tmem132d

Note: This product is for in vitro research use only