

# Taste receptor type 2 member 38

Catalog No: #AP88029



Package Size: #AP88029-1 50ug #AP88029-2 100ug #AP88029-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Taste receptor type 2 member 38
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	T2R38,T2R138,mT2R31,T2r31,Tas2r138
Accession No.	Uniprot ID:Q7TQA6
Uniprot	Q7TQA6
GeneID	387513;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Tas2r38

Note: This product is for in vitro research use only