

Translational activator of cytochrome c oxidase 1

Catalog No: #AP88030

Package Size: #AP88030-1 50ug #AP88030-2 100ug #AP88030-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Translational activator of cytochrome c oxidase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Coiled-coil domain-containing protein 44, Translational activator of mitochondrially-encoded cytochrome c oxidase I, CCDC44, PRO0477
Accession No.	Uniprot ID: Q9BSH4
Uniprot	Q9BSH4
GeneID	51204;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TACO1

Note: This product is for in vitro research use only