

Kit ligand

Catalog No: #AP88035



Package Size: #AP88035-1 50ug #AP88035-2 100ug #AP88035-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Kit ligand
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Hematopoietic growth factor KL,Mast cell growth factor,MGF,Stem cell factor,SCF,c-Kit ligand,Kitl,Mgf
Accession No.	Uniprot ID:P21581
Uniprot	P21581
GeneID	60427;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Kitlg

Note: This product is for in vitro research use only