

# Semaphorin-6B

Catalog No: #AP88042



Package Size: #AP88042-1 50ug #AP88042-2 100ug #AP88042-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Semaphorin-6B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Semaphorin-Z,Sema Z
Accession No.	Uniprot ID:O70141
Uniprot	O70141
GeneID	84609;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Sema6b

Note: This product is for in vitro research use only