

SUMO-activating enzyme subunit 1

Catalog No: #AP88044



Package Size: #AP88044-1 50ug #AP88044-2 100ug #AP88044-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SUMO-activating enzyme subunit 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Ubiquitin-like 1-activating enzyme E1A,AOS1,SUA1,UBLE1A
Accession No.	Uniprot ID:Q9UBE0
Uniprot	Q9UBE0
GeneID	10055;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SAE1

Note: This product is for in vitro research use only