

## Serpín B6

Catalog No: #AP88046



Package Size: #AP88046-1 50ug #AP88046-2 100ug #AP88046-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Serpín B6   |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE  |
| Species Reactivity | Ms  |
| Other Names        | Placental thrombin inhibitor, Proteinase inhibitor 6, PI-6, Serpinb6a, Spi3 |
| Accession No.      | Uniprot ID:Q60854   |
| Uniprot            | Q60854  |
| GeneID             | 20719;  |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                      |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Serpinb6

Note: This product is for in vitro research use only