

Tryptophan--tRNA ligase, cytoplasmic

Catalog No: #AP88047



Package Size: #AP88047-1 50ug #AP88047-2 100ug #AP88047-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tryptophan--tRNA ligase, cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Interferon-induced protein 53,IFP53,Tryptophanyl-tRNA synthetase,TrpRS,IFI53,WRS
Accession No.	Uniprot ID:P23381
Uniprot	P23381
GeneID	7453;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:WARS

Note: This product is for in vitro research use only