

Probable helicase senataxin

Catalog No: #AP88048

Package Size: #AP88048-1 50ug #AP88048-2 100ug #AP88048-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Probable helicase senataxin |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Amyotrophic lateral sclerosis 4 protein, SEN1 homolog, Senataxin, ALS4, KIAA0625, SCAR1 |
| Accession No. | Uniprot ID:Q7Z333 |
| Uniprot | Q7Z333 |
| GeneID | 23064; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SETX

Note: This product is for in vitro research use only