

# Tyrosine-protein phosphatase non-receptor type substrate 1

Catalog No: #AP88053

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP88053-1 50ug #AP88053-2 100ug #AP88053-3 1mg

## Description

Product Name	Tyrosine-protein phosphatase non-receptor type substrate 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	SHP substrate 1,CD172 antigen-like family member A,Inhibitory receptor SHPS-1,MyD-1 antigen,Signal-regulatory protein alpha-1,Sirp-alpha-1,CD172a,MYD1,PTPNS1,SHPS1,SIRP
Accession No.	Uniprot ID:O46631
Uniprot	O46631
GeneID	327666;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SIRPA

Note: This product is for in vitro research use only